GE Healthcare Life Sciences IN Cell Analyzer 6500HS Standard Operating Procedure

A. Overview

IN Cell Analyzer 6500HS is specialized in performing repetitive imaging tasks on a microplate format. Data generated by the instrument is grouped into a folder containing TIFF images and an ".xdce" file directory storing positional metadata of the TIFF image acquired. The user can choose to view the images one by one but if one wish to analyse the images together, please keep tiff and .xdce files together.

Before an imaging session. Please bring along the following:

- 1.1. Dimension chart of the plate the cells were seeded onto.
- 1.2. You have a map of the microplate you wish to image. Including positive and negative labelling control.
- 1.3. Excitation and emission Spectra of the dyes that have been used to label the cells if possible. Otherwise, name of the dyes. You will need a nuclear stain in order to count cell with automated procedure.
- 1.4. Portable data storage device of >100 GB of space.

B. Turn ON procedure

Turn on button (1) on machine beside gas inlet tubes. Turn on computer (2).





Wait until orange and red indicator light are extinguished.

Activate IN Cell Analyzer 6500 v7.3 software on desktop.



ell Analyzer 6500HS



Live cell incubator (**optional**) if your experiment require live cell incubation, please follow these steps.

Prerequisite: IN Cell Analyzer 6500 software is activated. (Critical*) Retrieve EC cover box from cabinet. Put into a BSC hood. Open lid and invert the EC cover turn on UV

irradiation for least 15 min.

Open valve on gas cylinders.

(Critical*) Turn both second regulator for compressed dry air (CDA) and CO₂ to 25 psi (black marking). **If your experiment require hypoxic incubation. Contact corefac@hku.hk for advanced setup at least one week before experiment.

[Optional] Turn N₂ gas supply regulator to 25 psi.

Bring plate into BSC hood remove plate cover and change into EC cover. **Keep you original cover in container to prevent contamination.

Connect first gas outlet then EC cover power cable. Activate incubation control software by clicking on the Environmental control icon on software.

| T. SAL |
|--------|

| | 3 | • | |
|-------|------------------------|-------------------------------------|-----------|
| •.(| Results • J≡ • | J | |
| | Environmo New style | ental Control controller with ga | s mixe |
| indow | | | - 0 |
| 2 26 | 5°C 🏼 | 47.1°C | \$ |
|) 15 | 0% | | |

| A window labelled "Env control window" should appear. | PEnv control window |
|---|-----------------------------------|
| | |
| | |
| 0.0 | |
| | Flow rate (I |
| \$O | Running (Heater Lid disconnected) |
| .0, | |
| | |
| | |
| \mathbf{C} | |
| | |
| | |

| | 26.5°C 🛛 47.1° | °C 🛸 |
|---------------------|----------------------------|---------|
| 0N 02 OFF | 15.0% | 1 |
| | °% | Ø |
| | Flow rate (humidified gas) | Ø |
| | | ъ 🖰 🔹 🤅 |

C. IN Cell Analyzer interface



<u>Top Left</u> panel for functions.

<u>Right</u> panel for image preview. Scroll to zoom.

sonthe solution

Lower left for plate map: scroll to zoom and double click a spot WITHIN plate area to move the lens to that position.

D. Create acquisition protocol

Prerequisite:

You have a list of dyes

You have information about the plate (minimum: plate product number; optimum: a manufacturer info sheet) Click start a new protocol in the software interface.





D2__Assign a name for the protocol. Name order [PI initial's]_[Your name]_[10 character description].

D3__Done. Please be reminded to click "save" after you have made major changes to the protocol or before starting



🔊 🇳 I

D4_ to run a protocol, Press "Scan" and then assign a folder to store your images. Please note that a folder will be created with the protocol's name and then image folder as sub-directory.

Acquire Preview Continuous Scan Stop

***Chapters E,F,G,H,I describes modules you can toggle during protocol set up (J,K ARE OPTIONAL).

E. Verify Plate dimension

E1__Select "Dashboard" page if you are not already on that page. Most of the setting will be done on this page.



E2__Click "Eject" to open the machine slide door for loading a plate.

E3__Insert microplate with A1 on top left corner of the sample frame.

E4___Click "Load" to close slide door.



E5__Choose correct plate from list. ***Important*** If you cannot find your plate catalogue number please contact FCF technical staff for assistance.

E6_ Choose objective you wish to use. Please see below guide for reference.

- 4X cell cluster counting
- 10X individual cell counting (large amount)
- 20X single cell morphology quantification
- 40X sub-cellular morphology quantification

E7__Verify plate bottom height with the lens you are going to image with. Verify again if you now decided to use a different objective to image.

E8____If the verify window indicates "100% of max" then decrease laser power by 4% and re-test. If blue line correctly identifies two peaks then click "Apply measured parameters".





| Plate/Slide: | Nunc 1251752 384-well plate | • | | Verify | |
|--------------|---------------------------------|---|---|--------|---|
| Obiective: | Nikon 4X/0.20. Plan Abo. CFI/60 | - | 6 | SAC , | , |

Objective: Nikon 4X/0.20, Plan Apo, CFI/60

- (4.4.4

•

E9 decrease laser power if needed in step E8 by

adjusting slider position and then click on pencil button

Power Level Power: 2% Power: 2% Power: 2% 1% 1% Use Channel Offsets

Another "Brightfield" channel available.

F. Channel selection

F1___ Click on dashboard. Plus sign. F2___ Select the color you wish to image.

Fluorescent channels available:

if required)

"ORANGE"

"FarRed"

"BLUE" "GREEN"

F3__ Click on camera icon check intensity. Adjust "Exposure" and "Laser Power" Accordingly. Please check fluorescent intensity below preview window. Maximum should be around 10,000 to 20,000 intensity units for fixed sample. Decrease laser power so that you image intensity unit of ~7,000 – 9,000 for live cell timelapse imaging to decrease phototoxicity. Quickly decrease laser power if maximum intensity reaches 65,535 units.

Min: 438 | Max: 47398 | Mean: 27925 | %Optimum: 90% | Z: 2934 Focus Finder

F4_____if you require confocal mode for deep depth of focus opticall

sectioning. Contact technical staff for aperture calibration now. Aperture calibration is done every morning.

(Please hover your cursor on top of the channel to read filter specifications

F5___if IRIS confocal mode:close down aperture to 1 AU to 2 AU (AU = airy units).

F6___ if EDGE confocal mode: select "Edge confocal" and close down aperture to 1AU to 2AU.

F7___monolayer cell culture works the best with "open aperture" selected. [N.B.] settings for F5 to F7 can be found here:

Fluorescent image modes: "2D", "3D" and "max intensity projection" Brightfield image modes: "2D", "DIC", "PH" (Phase contrast) and "3D" (please activate "3D" if you are doing IN Carta V1.10 volumetric analysis for brightfield, if activated.)

F8_____if your specimen has multi-colour labels, please find a site with all of the colours and click "Auto-offset".**

G. Field placement

| Dright leid Sett | rings | | | |
|------------------|----------------------|-----------|--|--|
| Expos | ure: 0.0 | 02000 🚖 | | |
| AF Off | iset: | 0.00 | | |
| Image Me | ode: 2-D | • | | |
| | Fran | ne Avg. 1 | | |
| 🗹 Open Aperl | ture 📄 Edge Confocal | | | |
| Laser Pov | wer: | 096 | | |
| Aperture (| AU): | 3.13 🔹 | | |
| I AF | | | | |
| Auto Offset | | | | |







| A | |
|---|---------------|
| 8 | |
| c | |
| D | |
| ε | 0000000000000 |
| F | 000000000000 |
| G | 000000000000 |
| н | 000000000000 |
| | |

Click preview on top-right corner.

| Acquire | Preview | Continuous | Scan | Stop |
|---------|---------|------------|------|------|
| | | | | |

| Imaging | Point Selection |
|---------|-----------------|
|---------|-----------------|

H. Focus strategies

| Click preview on top-right corner. |
|--|
| |
| Click "Point selection" to define scan areas. Imaging Point Selection |
| H. <u>Focus strategies</u> |
| H1_You can change focus position manually with "Focus Finder" tool. Click on Free under manual bar. You will be able to find Focus finder window from the |
| Focus Finder |
| |
| -1480 Protocol Z _o ; <u>3148</u> 2934 5020 |
| H2_You can activate autofocus once by clicking on each focus type icons: |
| H3_Add autofocus strategy to sampling rule by check boxes. You can test the autofocus at different positions on the plate with "AF" button are. |
| ✓ Software Autofocus ● First channel only ● Every channel On every 1 → field(s) of every 1 → well(s) |
| H4_ General guideline on which focus strategies to use: Adherent cell, use laser autofocus + auto offset. Z-stack with fixed positions, deactivate both autofocus and use "initial focus" define position Unpredictable heights (e.g. spheroids in 3D matrix), use software autofocus. |

I. Data transfer policies

I1_ After your protocol has finished running. There should be a folder containing: 1) an .xdce file; 2) tiff images; 3) a folder named "thumb". Be sure to transfer the whole folder.

I2_Go to computer (file explorer) select "IN-CARTA" under "Network"

 $\ensuremath{\mathsf{I3}}\xspace$ transfer your folder to one of the IN_Carta drives.

Folder naming rule: IN_Carta drive# / [Your lab PI's initials + "Lab"] / [Your name] / [Your image folder]

J. Z-stack setup (optional)

This function would be the most useful if used together with EDGE or IRIS confocal fluorescent imaging mode. Please refer to last section. Make sure you have selected "3D" for every channel you wish to acquire z-stack.

J1___On left menu area, select "z-stack setup"

J2___ use "FF" focus finder tool to focus onto the structure you wish to capture.

J3___ decide on z-stack optical section interval. The software will provide "Depth-of-field" in μ m for each objective. Maximum z-stack interval available is 50 μ m.

J4__select current z-position _____ to set as 3D focus result height (mid-point of z-stack).

J5____ use up and down arrow to define your "Start slice" and "End slice" positions. [N.B.] The z-range is governed by z-interval and number of slices. If z-interval is changed, slice number will remain constant.

| End Slice: | |
|-------------|--|
| Start Slice | |

F_F

1

20.70 🚔

(Depth-of-Field: 20.6 µm)

s 🧭 関

Dashboard
 Protocol Info
 Plate/Slide
 Objective Lens
 Fields
 Slide Imaging
 Channel Settings
 Focus Options
 Z-Stack Setup
 Image Processing
 Plate Heaters
 Liquid Handling
 Time Series
 Acquisition Options

 F_F

Z Step:

-entre (0)

| K. <u>Time series setup (optional)</u> | |
|---|--|
| K1_Select "Time Series" on menu. | Focus Options Z-Stack Setup Image Processing Plate Heaters Liquid Handling Time Series Acquisition Options Smart Scan |
| K2_check acquire time series button K3_determine "Time interval". It might be easier if "Display hh:mm function is turned on Display hh:mm:ss.ss | :ss.ss" Time Interval: 1:00 |
| K4_ Determine the number of time points required. | Number of Time Points: 4 |
| K5_If more than one well needs to be imaged. Please switch Mode "Multi-well". K6_click "Add" | to |
| K2 🗹 Acquire time series Starting Time Point: | 0 🔹 Estimate Burst Duration |
| K3 ☑ Display hh:mm:ss.ss K3 Time Interval: | 1:00 |
| K4 Number of Time Points: | 4 Add K6 Delete |
| ✓ Incubate between time points ✓ Refocus at each Time Point | |
| centre | |

L. <u>Turn OFF procedure</u>

L1_Live cell incubator shutdown (Optional, if you are working with fixed samples skip forward to step "L7")

L2_Turn off secondary regulator of compressed dry air and CO2. Keep airflow at 100 \mbox{cm}^3 / minute.

L3_(Critical**) Wait until both gas pressure decrease to 0.

L4_Turn off all functions in "Env control window". Then close Env control.

L5_Press "Eject" button in IN Cell software

L6_Unplug gas tubing and heating connector. Return the gas tubing to its holder.

L7_Turn off IN Cell Analyzer.

***Please do this step AFTER live cell incubator has been shut down (if applicable).

- L8_ Eject any plate that is still left on sample holder.
- L9_Go to: Application-> Hardware-> Shutdown Instrument

sonthe top another and the sonther and the son

- L10_Transfer data you need from IN Carta offline station.
- L11_Turn off computer.

clences