

# COY In-vitro Hypoxia system

<b>Specifications</b>	
<b>Cell Chamber:</b>	Two Chambers with controller to monitor and control O <sub>2</sub> level (0.1%-21%) Chamber Dimension – 370mm x 410mm x 220mm
<b>Incubator:</b>	NuAir cell incubator set at 37° C
<b>Gases supply:</b>	5% CO <sub>2</sub> in N <sub>2</sub> and 5% CO <sub>2</sub> in O <sub>2</sub> compressed gas
<b>Application</b>	
Chronic hypoxic tissue culture	